


# SOZO® Architecture Diagram



## TECHNICAL SPECIFICATIONS

- Bluetooth v 2.1 encryption standard based communication of the Tablet with the SOZO Device
- Internet enabled Wi-Fi for the Tablet to communicate with ImpediMed servers on Port 443

## SECURITY AND PRIVACY OF DATA

To provide safe and reliable access to the data, SOZO Digital Health Platform offers the following security and privacy features.

- Requires SOZO users to have complex application passwords
- Ability for the end users (clinician and administrator) to reset their own password
- Configurable password lockout period and password expiration time
- Ability to create multiple administrator accounts
- No data resides on the SOZO Device or the Tablet
- Data in motion and Data at rest are encrypted
- Availability of detailed user-logs via user interface
- High availability and enhanced encrypted data back-up
- Segregated virtual private network for different objects in the architecture
- Access to ImpediMed SOZO infrastructure restricted to authorised users with 2-factor authentication
- Have implemented the following security and privacy controls
  - a. Security Standards
  - b. Administrative Safeguards
  - c. Physical Safeguards
  - d. Technical Safeguards
  - e. Organisational Requirements
  - f. Policies and Procedures
  - g. Notification to the Secretary
  - h. General Rules; Uses and Disclosures of PHI
  - i. Organisational Requirements; Uses and Disclosures

## impedimed®

ImpediMed Limited  
Unit 1, 50 Parker Court, Pinkenba, QLD 4008 Australia  
Phone: +61 7 3860 3700  
Toll Free: 1800 638 477  
Email: [enquiries@impedimed.com](mailto:enquiries@impedimed.com)  
[www.impedimed.com](http://www.impedimed.com)

©2019 ImpediMed Limited. ImpediMed and SOZO are registered trademarks of ImpediMed Limited.